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
Amerikansk patentnr. 5,209,754

VCARE er designet af Dr. Prabhat K. Ahluwalia, Little Falls, NY.

**VCARE®****Vaginal-Cervical Ahluwalia's Retractor-Elevator**

Reference # 60-6085-200A, -201A, -202A, -203A

OD Size: 5 mm (0.2 in)

Length: 42 cm (16.4 in) including Handle

Disposable Single-Use


Read and become familiar with all instructions, warnings, and cautions in the Package Insert before using this product.

Sterile: This product is sterile in the unopened, undamaged package.

Note: This product is not made with natural rubber latex.

Description:

VCARE is a disposable, single-use device for manipulation of the uterus and cervix in surgical and diagnostic procedures. VCARE consists of an insulated manipulator tube having an inflatable balloon at its proximal [to patient] end and an anatomically configured cannula/handle for maintaining proper attitude of the uterus at the distal end. (Figure #1) The insulated manipulator tube is marked with reference graduations from the proximal end (centimeters). The graduations are provided as a guide for comparison to a graduated uterine sound. VCARE incorporates a system of cup-like elevators to provide manipulation of the uterus, and retraction and elevation of the cervix. There are four (4) size variations of cervical cups - Small (S), Medium (M), Large (L), and Extra-Large (XL). (Figure #2) The cervical cups are green and provide a colpotomy guide. Sites for suturing are located on the cervical cup.

Figure #1**Indications:**

The CONMED VCARE Retractor/Elevator is indicated for manipulation of the uterus, and injection of fluids or gases during laparoscopic procedures such as Laparoscopic Assisted Vaginal Hysterectomy (LAHV), Total Laparoscopic Hysterectomy (TLH), minilap, laparoscopic tubal occlusion or diagnostic laparoscopy and also maintains pneumoperitoneum by sealing the vagina once colpotomy is performed.

Contraindications:

This device should not be used when insertion of intrauterine retracting devices is contraindicated.

This device is not intended for use as a Laser Backstop.

The use of an intrauterine device may be contraindicated for cervical or uterine pathology. VCARE should not be used on patients who are, or are suspected of being pregnant.

Directions for Use:

Read and become familiar with all instructions and cautions before using this device.

A. VCARE SIZE SELECTION:

- Visually examine the patient's cervix and assess its volume.
- Choose a VCARE which has a cervical cup that most closely approximates size of cervix. (Figure #2)

Figure #2

Reference Dimensions			
VCARE Size Guide	Cervical Cup Volume	Cervical Cup Diameter	Cervical Cup Designation
Small, REF 60-6085-200A	9.7 cm ³ (0.59 in ³)	32 mm (1.26 in)	S
Medium, REF 60-6085-201A	14.3 cm ³ (0.87 in ³)	34 mm (1.34 in)	M
Large, REF 60-6085-202A	20.7 cm ³ (1.26 in ³)	37 mm (1.46 in)	L
Extra-Large, REF 60-6085-203A	26.7 cm ³ (1.63 in ³)	40 mm (1.57 in)	XL


B. INSERTION OF VCARE INTO THE UTERINE CAVITY:

- Remove VCARE from its sterile packaging and inspect for damage caused by shipping. Discard if any damage is noted. Draw 7-10 cc of air into a standard syringe and insert it into the Luer connection at the end of the pilot balloon. Test the intrauterine balloon by injecting air with the syringe and checking to see if balloon remains inflated. If balloon does not remain inflated, DO NOT USE. Discard and select another VCARE unit. After the successful balloon test, deflate the balloon by removing all air with the syringe.
- The patient should be placed in the desired surgical or diagnostic (lithotomy) position.
- Following surgical technique, and ensuring good visualization of the cervix, the surgeon should grasp the anterior cervical lip with an Allis clamp or similar (toothed) instrument.

- The direction of the uterine canal and the depth of the uterine cavity are determined by use of a blunt probe or graduated uterine sound. The graduations are provided as a guide for comparison to a graduated uterine sound.
- If necessary, dilate the cervix to accommodate the 5 mm (0.2 in) diameter of the manipulator tube.
- If suturing of the cervical cup is desired, sutures should be placed prior to device insertion.
- Lubricate intrauterine balloon and tip with a sterile surgical lubricant.
- Carefully insert the proximal tip of VCARE through the Cervical Os until the balloon is in the desired position within the uterine cavity. NOTE: If the uterus is retroverted, rotate the VCARE 180° prior to inserting the intrauterine balloon. Once inserted, rotate the VCARE 180° to antevert the uterus.
- Remove the clamp from the cervix.
- Stabilize the distal (external) end of VCARE to prevent accidental dislodging or forward movement of the device until fully positioned.
- Reattach the syringe to the Luer adapter on the pilot balloon and inflate the intrauterine balloon with at least 7 cc of air. Maximum inflation is 10 cc of air. (Figure #3)
- Remove syringe to prevent spontaneous deflation of the balloon by backward pressure.

POSITIONING THE CERVICAL/VAGINAL CUPS:


- Slide the cervical cup along the manipulator tube until its outer edges surround the vaginal fornix. Ensure proper positioning of this cup before proceeding. NOTE: If suturing the cervical cup in place, thread suture through suture slots prior to cervical cup positioning. Once cervical cup is in place, suture can be tied off.
- A. Slide the vaginal cup along the manipulator tube until it meets the back edge of the cervical cup and is completely seated against the vaginal walls to ensure maintenance of pneumoperitoneum as required. B. Slide lock down to vaginal cup and secure in place by turning thumbscrew clockwise until tight. (Figure #3)
- To ensure that the intrauterine balloon has not ruptured during insertion, check the tautness of the pilot balloon. (See WARNINGS)
- Apply gentle traction to the manipulator tube to check that VCARE is secure and that the uterus is grasped.

Figure #3


+1 (315) 797-8375 • FAX +1-800-438-3051
Kundeservice (USA) +1-800-448-6506
E-mail: info@conmed.com

© CONMED Corporation, trykt i USA

Kontakt CONMED International Sales ved internationale order eller forespørgsler
+1 (315) 797-8375 • FAX +1 (315) 735-6235
www.conmed.com

**C. REMOVAL AND DISPOSAL OF VCARE:**

- Re-attach the syringe to the Luer connector at the end of the pilot balloon; fully aspirate the air from the intrauterine balloon to deflate. This will allow the intrauterine balloon to be removed from the uterus.
- Unlock the locking mechanism by turning the thumbscrew counter-clockwise (anti-clockwise) and retract to the handle.
- A. Swipe finger around the edge of the vaginal cup to separate the tissue from the cup to prevent tissue damage. B. Fully retract the vaginal cup to the handle. (Figure #4)
- Carefully remove the device from the vagina. DO NOT use excessive force to avoid traumatizing the vaginal canal.
- Upon removing VCARE, the surgeon should visually inspect the VCARE device, and the patient, to make sure that the entire VCARE device was properly removed and that no components or fragments of these components were retained in the patient. There are 5 parts/components to the VCARE Cervical Elevator Retractor. These are: 1) the balloon; 2) the forward "cervical" cup; 3) the back or vaginal cup; 4) the locking assembly with thumb screw; 5) the metal shaft and handle with balloon inflation valve.

Figure #4**INTRAUTERINE INJECTIONS:**

VCARE has a channel for injection of appropriate dyes or marking materials. Remove the cap and insert the syringe on the Luer connector at the rear of the handle. Use slow and steady pressure for injection of prescribed fluids. (See WARNINGS)

WARNINGS:

- For safe removal of the VCARE device from the patient, follow instructions in Section C. Failure to separate the vaginal cup from the tissue may result in detachment of the cervical cup and/or patient injury. (Figure #4)
- Prior to device removal, ensure the locking mechanism is released via the thumbscrew and swipe a finger around the edge of the vaginal cup to separate the tissue from the cup to prevent tissue damage. (Figure #4) DO NOT use excessive force upon device removal to avoid traumatizing the vaginal canal and/or component detachment.
- Vaginal delivery of a large uterus may result in patient injury. Methods should be used to reduce the size of the uterus prior to removal through the vaginal canal.
- Upon removing VCARE, the surgeon should visually inspect the VCARE device, and the patient, to make sure that the entire VCARE device was properly removed and that no components or fragments of these components were retained in the patient. There are 5 parts/components to the VCARE Cervical Elevator Retractor. These are: 1) the balloon; 2) the forward "cervical" cup; 3) the back or vaginal cup; 4) the locking assembly with thumb screw; 5) the metal shaft and handle with balloon inflation valve.
- This device should be used only by surgeons trained in proper techniques for uterine surgery, diagnostic procedures, gynecological pelvic anatomy, and placement of intrauterine retracting instruments. Read and become familiar with all instructions, warnings and precautions in this Package Insert before using this device. As with all intrauterine devices, uterine perforation and subsequent bleeding is a possibility. The surgeon should examine the patient carefully for any indications of uterine bleeding and take appropriate clinical steps.
- Improper use of this or any intrauterine instrument can result in perforation of the uterine wall and subsequent bleeding.
- NEVER use this device with the intrauterine balloon deflated.
- To avoid potential injury to the uterine wall or expulsion of the device from the uterus do not underinflate the intrauterine balloon. Inflation with 7 to 10 cc of air is recommended to compensate for injected air trapped in the dead space of the pilot balloon and inflation tube. DO NOT exceed 10 cc.
- DO NOT use any liquids to inflate the intrauterine balloon. Use of fluid could potentially over-distend and burst the balloon since fluid is not compressible.
- DO NOT inject dyes or marker fluids into the pilot balloon. Dyes or marker fluids should only be injected through the Injector Port on the rear of the handle. Do not attach other devices or lines to the pilot balloon.

- DO NOT inject dyes/marketing fluids rapidly as high intrauterine pressures can develop. Injection of such fluids can result in expulsion of the device from the uterus, spasm of the Fallopian Tubes or vascular extravasation. Use slow and steady pressure for injection of any appropriate dye or marker fluid.

- DO NOT attach other devices or lines other than the syringe to inject contrast medium to the injector port.
- VCARE must be inserted into the uterus along the direction of the uterine cavity, i.e. anteriorly in anteverted uteri, posteriorly in retroverted uteri, to reduce potential for uterine perforation. Perform pelvic examination to determine the direction of the uterus and then sound the cavity to assess the uterine depth. Extra caution should be exercised in the case of a very small uterus. Do not use VCARE if the uterus sounds to less than 4.0 cm (1.6 in).

- DO NOT use VCARE as a uterine sound.
- VCARE is not recommended for use in a large, post-partum uterus. Physicians must exercise sound judgement and care, should they choose to use this device. Movements of the manipulator within the uterus may cause lacerations within the uterine wall and bleeding.

- Check the Pilot Balloon frequently to ensure inflation of the intrauterine balloon. If the intrauterine balloon ruptures, the pilot balloon will not feel firm when squeezed between your fingers. If the intrauterine balloon has ruptured, stop all manipulation immediately, remove VCARE and replace it with a new VCARE unit. Use of a ruptured balloon eliminates the protection to the uterine wall and may increase the potential for uterine perforation.
- The cervical/vaginal cup locking mechanism must be locked at all times when using VCARE for uterine manipulation. If lock inadvertently becomes open, secure it immediately before proceeding.
- In laparoscopic supracervical hysterectomy, avoid contact between VCARE and dissecting instruments during uterine dissection and excision to avoid risk of patient injury. Damage to VCARE, including rupture of the intrauterine balloon, may occur.
- If bending VCARE is necessary for use with robotic instruments, then all device systems must be tested before use to ensure functionality.

- Following use of this device and where the uterus has not been removed, the surgeon should check for any indications of uterine perforation and subsequent bleeding. Should uterine perforation and/or bleeding occur, immediately follow the clinically indicated procedures to establish hemostasis.
- DO NOT reuse this device. The ability to effectively clean and re-sterilize this single use device has not been established and subsequent re-use may adversely affect the performance, safety, and/or sterility of the device.

- This product contains phthalates that may leach out into the patient during use. Animal studies have shown that high-level exposure to DEHP may affect male fertility and reproductive development. Care should be taken to minimize exposure to children and nursing women.

CAUTIONS:

- Test the intrauterine balloon for proper inflation and function before insertion; do not insert VCARE with intrauterine balloon inflated.
- Lubricate the intrauterine balloon and end of the tube before insertion.
- Check for the need to dilate cervix before the insertion of VCARE. The Cervical Os should be a #13-14 Hank size or at least 5 mm (0.2 in) in diameter before insertion to allow for free passage and to prevent injury to the cervix, potential for uterine perforation, or risk of tearing the intrauterine balloon.
- As with all intrauterine devices, uterine perforation and subsequent bleeding is a possibility. The surgeon should examine the patient carefully for any indications of uterine bleeding and take appropriate clinical steps.
- When injecting any liquid media, closely follow the manufacturer's directions for use accompanying the product.
- Use of VCARE in hysterosalpingography and Rubin's Test should be completed under anesthesia.

U.S. Patent # 5,209,754

VCARE was designed by Dr. Prabhat K. Ahluwalia, Little Falls, NY.

For International orders or inquiries,
Please contact CONMED International Sales
+1 (315) 797-8375 • FAX 1-800-438-3051
Customer Service (USA) 1-800-448-6506
e-mail: info@conmed.com

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© CONMED Corporation, gedrukt in de VS