# PlumePen® Ultra

## Surgical Smoke Evacuation Pencil

**Intended Use:** PlumePen® Ultra is designed for general electrosurgical applications, including cutting and coagulation, and for removing smoke generated by electrosurgery when used in conjunction with an effective smoke evacuation system. The pencil enables the operator to remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the operative site for the desired surgical effect to remove tissues and control bleeding by means of high-frequency electrical current.

**Indications for Use:**

- a. To remove smoke plume from the surgical site.
- b. To remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the target tissue for the desired surgical effect.


**Figure 1.** Open product packaging and remove PlumePen® Ultra assembly. Connect the tubing provided with PlumePen® Ultra to a surgical smoke evacuator or other appropriate vacuum source; example shown. (Figure 1). For connection and operation to a specific surgical smoke evacuator, please refer to the specific instructions for use that accompany that system for appropriate set-up and use.

**Figure 2.** Plug the 3-prong power cord into the electrosurgical generator of your choice, (Figure 2). Confirm that all power settings on the generator are appropriate for the procedure being performed.


The electrosurgical generator's intensity should be set as low as is necessary to achieve the desired effect. PlumePen® Ultra is a monopolar electrode, the use of a dispersive electrode is required to prevent burns/injury to patient.

Please refer to electrosurgical generator user manual and dispersive electrode instructions for use for additional instructions.

**Figure 3.**  
**3a.** The PlumePen® Ultra is equipped with several key features for the surgeon's convenience. Surrounding the electrode blade is a translucent plastic tube that allows the blade and tube to be positioned at varying lengths to most effectively capture surgical smoke plume as it is created.  
**3b.** Rotate the CAM locking mechanism to the unlock position.  
**3c.** Extend Capture Port / Blade to desired length. Please note Capture Port and Blade move together.  
**3d.** Turn the CAM locking mechanism back to the lock position and ensure that the tube is securely locked in place.

Avoid contact with blade and buttons when adjusting the capture port. Keep the active electrodes clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning. Injury to operating room personnel may result.


**Figure 4.** The PlumePen® Ultra has CUT and COAG buttons located on the top of the electrosurgical pencil body. The YELLOW button is for cutting and the BLUE button is for coagulation, (Figure 4).

After use, turn off smoke evacuation unit and electrosurgical generator. Discard the entire PlumePen® Ultra assembly (Pencil, Blade, Tubing, and Power Cord) per facility protocol for contaminated items.

**Figure 5.** To optimize the life of your smoke evacuation filter, it is recommended that a Buffalo Filter Surgical Smoke Plume Evacuator be used in conjunction with an EZLink® Automatic Activation Device (Model EZLink01), (Figure 5).

The EZLink® will automatically turn the Buffalo Filter Surgical Smoke Evacuator on or turn it off when either the CUT or COAG function is activated/de-activated on the PlumePen® Ultra. Contact CONMED for more details on the EZLink®.

**Figure 6.** If necessary to remove blade, **6a.** Unplug the pencil **6b.** Ensure that CAM is in the lock position **6c.** Use a surgical clamp and pull blade forward **6d.** Replace with blade of choice **6e.** Confirm that blade is completely inserted and secure before activating pencil. Never force the blade into the pencil.

**CAUTION:** We do not recommend re-using the original blade after it has been removed.

**CAUTION:** Intended for use with a maximum voltage of 5 kVpk.

**DO NOT USE** in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.

Sterile unless packaging is damaged or any seal is broken. Do not use electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids, or objects, or in the presence of oxidizing agents, as fire may result. When not in use, the active electrode should be placed in a clean, dry, non-conductive safety holster. Inadvertent contact with the patient may result in burns. Contact with drapes or linens may cause a fire.

DO NOT activate the instrument when not in contact with target tissue, as this may cause injuries due to capacitive coupling.

PlumePen® Ultra is designed and intended only to be used with an electrosurgical generator that has been tested to the IEC 60601 standard. Please refer to generator to ensure compatibility. PlumePen® Ultra is not a fluid removal device, therefore, it should not be used for such application.

ASPIRATE fluid from the area before activating the instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient.

For procedures where visualization may be impaired, be alert of these potential hazards:

- The electrode tip may remain hot enough to cause burns after the current has been de-activated.
- Inadvertent activation or movement of the activated electrode outside of the field of vision may result in injury to the patient.
- Localized burns to the patient or physician may result from electrical currents carried through conductive objects. Current may be generated in conductive objects by direct contact with the active electrode, or by the active accessory being in close proximity to the conductive object.

For single use only. Do not reuse, reprocess, or re-sterilize. Reuse, reprocessing, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness, or death.

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