

- Step 3 Use soft bristle brush to manually clean the Battery with the cleaning solution.
- Step 4 Ensure areas containing crevices are scrubbed thoroughly.
- Step 5 Wipe off detergent thoroughly with lukewarm tap water.
- Step 6 Perform visual inspection to determine if debris is removed.
- Step 7 Repeat cleaning as necessary to obtain a visually clean Battery.

Chemical Disinfection

Disinfectants should be prepared and used according to the manufacturer’s recommendations. It is recommended that the chemical disinfectant be wiped off with tap water.

- 70% Isopropyl alcohol
- 10% Bleach (sodium hypochlorite solution)

Warnings and Precautions

- Examine the shipping carton and Instrument for signs of shipping damage. Note any shortages, breakage, or apparent damage, retain the evidence, notify Customer Service or Distributor immediately and replace with a new Instrument. Do not use a damaged product.
- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Instruments for minimally invasive procedure may vary in diameter from manufacturer to manufacturer. When such Instruments and accessories from different manufacturers are employed together in a procedure, verify their compatibility prior to procedure.
- Do not use the Instrument if the shaft is visibly bent.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- The Instrument must be disposed after procedure once the package is opened.
- The Instrument is designed, inspected, and manufactured for single procedure only. Do not reuse, reprocess or resterilize the Instrument as it may compromise the structural integrity of the Instrument, and/or lead to Instrument failure that in turn may result in patient injury, illness, or death.
- Reusing the Instrument may create risk of contamination, infection, or cross-infection, including, but not limited to, the transmission of infectious diseases, which may lead to injury, illness, or death.
- Do not load the Instrument more than 16 times. The Instrument can fire for a maximum number of 16 times. Use of staple line reinforcement material may reduce the maximum number of firings.
- After removing the Shipping Wedge, observe the surface of the Reload. The Reload Units must be replaced with another Reload Units if any staple tray is visible. (If staple tray is visible, the Reload may not contain staples.)
- Do not articulate when the jaws are closed.
- When selecting the Reload Units, careful consideration should be given to existing pathologic conditions as well as any pre-surgical treatment, such as radiotherapy, that the patient may have undergone. Certain conditions or preoperative treatments may cause change in tissue thickness that would exceed the indicated range of tissue thickness for the standard choice of Reload Units.
- Do not use hospital autoclaves to sterilize or disinfect Battery and the Instrument.
- Use of any other type of battery other than the battery supplied with the Instrument may result in increased EMISSIONS or decreased IMMUNITY of the Instrument.
- Avoid using the Instrument adjacent to or stacked with another equipment. If it is necessary to use the Instrument adjacent or stacked with another Instrument, pay attention, and notice any abnormalities.
- Do not modify the Instrument without authorization from the manufacturer.
- Use of accessories other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this Instrument and result in improper function.
- If the hemostasis of the staple line cannot be clearly observed, do not continue using this Instrument.

The Instrument must be used in a specified electromagnetic environment. For more information, refer to **Guidance and manufacturer's declaration for EMC**. Failure to follow these instructions may cause the Instrument to malfunction.

- The Instrument cannot be operated under oxygen enriched environment.
- A notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to Reach Surgical, Inc. through Reachquality@reachsurgical.com and the competent authority of the Member State in which the user and/or patient is established.

Technical Parameters

- Battery ratings: DC 12V,1550mAh.
- The Instrument is not resistant to fluids ingress and is classified per IEC 60601-1 as IPX0
- The Instrument needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

Storage Requirements

Temperature: 5°C ~ 35°C
Relative Humidity: 0 % ~ 70 %

Air Pressure: 500 hPa ~ 1060 hPa

Transport Requirements

Temperature: -10°C ~ 54°C
Relative Humidity: 0 % ~ 70 %

Air Pressure: 500 hPa ~ 1060 hPa

Operating Environment Requirements

Temperature: 10°C ~ 40°C
Relative Humidity: 30 % ~ 75 %

Air Pressure: 800 hPa ~ 1060 hPa

Software information

The Instrument is controlled by an embedded software program with the version of software V01. The software program is intended to detect the current. If the current is too high, the software will cut off the circuit to avoid motor damage.

Expiration Date

The Instrument is sterilized by Ethylene Oxide. The expiration date is labeled on the package. Do not use this Instrument beyond its expiration date.

How Supplied

This Instrument is supplied sterile for single patient use. Discard after use.

Guidance and manufacturer's declaration for EMC

WARNING: Use of the Instrument adjacent to or with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used not closer than 30 cm (12 inches) to any part of the Powered Articulating Staplers and Reload Units, including cables specified by the MANUFACTURER. Otherwise, the degradation of the performance of this equipment could result.

NOTE: The EMISSIONS characteristics of the Instrument make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), the Instrument might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Electromagnetic interference will not prevent the Instrument from being electrically closed and fired.

The Instrument must be used in a specified electromagnetic environment. Follow the directions of the following tables when using the Instrument.

Guidance and manufacturer's declaration – electromagnetic emission			
The Instrument and Reload Units are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:			
Emission test	Compliance	Electromagnetic environment - guidance	
RF emission CISPR 11	Group 1	The Powered Articulating Staplers and Reload Units use RF energy only for its internal function. Therefore, its RF emissions are low and there is little possibility of producing interference to nearby electronic equipment.	
RF emission CISPR 11	Class A	The Powered Articulating Staplers and Reload Units are suitable for use in professional healthcare facilities.	
Harmonic distortion IEC 61000-3-2	N/A		
Voltage fluctuations and flicker IEC 61000-3-3	N/A		
Guidance and manufacturer's declaration – Electromagnetic immunity			
The Powered Articulating Staplers and Reload Units are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV ±8 kV,±15 kV air	±8 kV contact ±2 kV, ±4 kV ±8 kV,±15 kV air	The floor should be wood, concrete or ceramic. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients / bursts IEC 61000-4-4	±2 kV 100 kHz repetition frequency	N/A	Battery powered and no signal line >3m
Surges IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground	N/A	Battery powered and no signal line >30m or going out to outdoor
Voltage dips IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315°	N/A	Battery powered
Voltage interruptions IEC 61000-4-11	70% UT; 25 cycles at 0° 0% UT; 250 cycle		
Rated power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz	30 A/m 50Hz	The power frequency magnetic field should have the characteristics for use in a typical place in a typical commercial or hospital environment.
Note: U _i refers to the AC voltage of the power supply before the test voltage is applied.			
Guidance and manufacturer's declaration –Electromagnetic immunity			
The Powered Articulating Staplers and Reload Units are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:			

Immunity test	IEC 60601 test level		Compliance level	Electromagnetic environment - guidance		
Conducted disturbances induced by RF fields IEC 61000-4-6 Radiated RF EM fields IEC 61000-4-3	3V 0.15 MHz ~ 80 MHz 6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1kHz 3 V/m 80 MHz ~ 2.7 GHz 80% AM at 1 kHz		N/A 3 V/m 80 MHz ~ 2.7 GHz 80% AM at 1 kHz			
Guidance and manufacturer's declaration –Electromagnetic immunity						
The Powered Articulating Staplers and Reload Units are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:						
Immunity to RF wireless communications equipment (IEC 61000-4-3)						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380—390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430—470	GMRS 460, FRS 460	FM ±5kHz deviation 1kHz sine	2	0.3	28
710	704—787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810						
870						
930	800—960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
1720	1700—1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25, UMTS	Pulse modulation 217Hz	2	0.3	28
1845						
1970						
2450						
5240	5100—5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						
Guidance and manufacturer's declaration –Proximity magnetic fields immunity						
The Powered Articulating Staplers and Reload Units are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:						
Immunity test	IEC 60601 test frequency		Modulation	Immunity test level(A/m)		
Proximity magnetic fields immunity IEC 61000-4-39:2017	134.2kHz		Pulse modulation, 2.1kHz	65		
	13.56MHz		Pulse modulation, 50kHz	7.5		
	EN Sterilization batch					
	EN Peel Here					
	EN HDPE recyclable					
	EN Recyclable					
	EN Type BF Applied Part					
	EN Refer to instruction manual					
	EN Authorized Representative in the European Community					
	EN Do not use if package is damaged.					
	EN Do not resterilize					
	EN Manufacturer					
	EN Date of manufacture					
	EN Serial number					
	EN Batch code					
	EN Use-by date					
	EN Fragile, handle with care					
	EN Keep dry					
	EN Keep away from sunlight					
	EN Up					
	EN Do not re-use					
	EN Caution					
	EN Catalogue number					
	EN Storage temperature limit					

	EN Storage humidity limitation	
	EN Storage atmospheric pressure limitation	
	EN Single sterile barrier system	
	EN Country of manufacture	
	EN Medical device	
	EN Unique device identifier	
	EN Sterilized by Ethylene Oxide.	
 www.int.reachsurgical.com/support 		EN Consult instructions for use or consult electronic instructions for use