

Preparation Bag with Memory Wire

SB836, SB857



Description

Unimax Preparation Bag – the preparation bag with memory wire is supplied with a cylindrical tube for easy introduction through a standard 10, 11 or 12 mm trocar.

Indications

The Unimax Preparation Bag – preparation bag with memory wire – is intended for use as a receptacle for the collection and extraction of tissue, organs and calculi during general and laparoscopic surgical procedures.

Contraindications

This device should not be used when endoscopic intervention is contraindicated.

Adverse effects

Keep the bag away from sharp instruments or edges. Rupture of the bag may result in spillage of fluid or tissue specimens.

Preparation

1. Remove the preparation bag from the packaging.
2. Insert the cylindrical tube with the preparation bag through the desired 10, 11 or 12 mm trocar.
3. To release the preparation bag from the tube, simply push the inner insert (1) forward and the bag unfolds automatically.
4. Place the desired specimen in the bag. Make sure that the entire specimen fits within the boundaries of the bag. When all desired specimens are in the bag, close the bag. To close the bag, grasp the green mark (2) on the memory wire with a grasper and pull the distal end of the memory wire (3) with a second grasper until the bag is completely closed.
5. Remove the specimen bag at the end of the laparoscopic procedure. To remove the bag, use a grasper to pull the top of the bag into the distal end of the trocar. While securing the bag with the grasper, withdraw the trocar, bag and grasper through the incision site until the leading edge of the bag can be visualised. Continue to remove the bag through the incision by hand and under observation.
6. When the endoscopic procedure is completed, discard the device in accordance with local regulations.

Precautions / Warnings

- This device is EO sterilised and for single use only. Do not reuse or reprocess, in order to avoid irreversible harm such as cross-infection or injury caused by device failure.
- Do not use if the packaging is open or damaged.
- Store the device in a dry, clean and safe place. Keep the device away from flammable and combustible chemicals or materials.
- Please follow the general guidelines for the disposal of medical waste to avoid possible cross-infection, environmental contamination or injury caused by sharp waste.