



# New Generation Endo Stapling System

P2G 





## TECHNOLOGY P2G

- Double Gripping Technology
- Lateral diffusion of the fluids
- Better tissue compression

The P2G technology guarantees better stapling delivering better hemostasis and safety

## BENEFITS

- Enhanced anvil delivers better compression & manipulation of the tissue
- Anvil and cartridge with progressive clamping technology for better clamping and less squeezing of the tissue
- Hemostasis & safety

## FEATURES

Double Gripping Technology



Drop-shaped staple groove



Progressive clamping technology



Enhanced anvil





**Panther**<sup>®</sup>  
Healthcare

Rounded tip for a safer  
and easier blunt dissection



Double Gripping Technology



Progressive Clamping technology.  
Anvil and cartridge with non-traumatic sides

Stronger  
clamping force

Enhanced  
I-Beam

Fixed Anvil  
with Rounded Tip

Coating



# NEW GENERATION ENDO STAPLING SYSTEM

Product Name	Product	Description	Color
UNIVERSAL Endo Linear Cutter	CEAB-30	Endo Linear Cutter (Standard Handle)	
	CEAB-45	Endo Linear Cutter (Short Handle)	
	CEAB-60	Endo Linear Cutter (Extra Long Handle)	

Endo Articulating Reloading	CADF-30T	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF-30N	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF-45T	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF-45N	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF-45S	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload	● Gold
	CADF-45D	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload	● Green
	CADF-45R	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload	● Black
	CADF-60T	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload	○ White
	CADF-60N	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload	● Blue
	CADF-60S	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload	● Gold
	CADF-60D	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload	● Green
CADF-60R	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload	● Black	