


# ConMed® GI Cytology Brushes



For Maximum Sampling of Cells  
for Diagnostic Testing

Applied Ingenuity™

**CONMED™**  
ENDOSCOPIC TECHNOLOGIES

# ConMed® GI Cytology Brushes

## BRISTLED BRUSH MATRIX

- Maximizes cell sample

## POSITIVE STOP IN HANDLE

- Ensures accurate brush position

## THUMB RING

- Provides "single handed" control

## MULTI-PURPOSE AND MULTI-SIZED

- Accommodates Adult and Pediatric Gastrosopes, Duodenoscope and Colonoscope

## DISPOSABLE AND STERILE

- Minimize risk of cross-contamination

000104

000108

000109

000110

000196

Straight Handle

Ring Handle

## ConMed® Sheathed Cytology Brushes

*Sterile, Disposable, Latex-Free*

Catalog Number	Description	Endo Scope	Brush Diameter	Length	Sheath O.D.	Unit
000104	White, Ring Handle	Multi-Purpose	1.7mm	160cm	1.7mm	20/box
000108	White, Ring Handle	Pediatric Gastroscope	3.0mm	160cm	1.7mm	20/box
000109	White, Ring Handle	Duodenoscope	3.0mm	200cm	1.8mm	20/box
000110	White, Ring Handle	Colonoscope	3.0mm	230cm	2.1mm	20/box
000196	White, Ring Handle	Adult Gastroscope	3.0mm	160cm	2.0mm	20/box

## Cytology Brushes

*Sterile, Disposable, Latex-Free*

Catalog Number	Description	Brush Diameter	Length	Sheath O.D.	Unit
151	Blue, Straight Handle	2.5mm	160cm	1.7mm	20/box
151R	Blue, Ring Handle	2.5mm	160cm	1.7mm	20/box
152	Blue, Straight Handle	2.5mm	200cm	1.8mm	20/box
152R	Blue, Ring Handle	2.5mm	200cm	1.8mm	20/box
163	Blue, Straight Handle	3.0mm	160cm	2.4mm	20/box
163R	Blue, Ring Handle	3.0mm	160cm	2.4mm	20/box
158	Blue, Straight Handle	3.0mm	230cm	2.4mm	20/box
158R	Blue, Ring Handle	3.0mm	230cm	2.4mm	20/box



## Biopsy Channel Accessory Introducer

*Non-Sterile*

Catalog Number	Description	Unit
000160	Accessory Introducer	6/box

The cannula/trocar device is used to open up the biopsy/valve that is sometimes restrictive and non-cooperative with tapered ERCP cannula, cytology brushes, manometry catheters, injection needles or polypectomy snares.

### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
Utica, NY 13502  
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com