



GENMED
ENTERPRISES

Laparoscopic

PRODUCT RANGE

Contents

- 1 Monopolar instruments
- 2 Disposable clip applicator
- 3 Veress needle
- 4 Retrieval bags
- 5 Suction irrigation

Monopolar Instruments



DETAILS

Full Range of Jaws Available

Finger Rotation Wheel Provides 360* Movement Colour Indicators

For Different Jaws

Fully Insulated Shaft

On/off Ratchet Mechanism to Securely Place Jaws

PRODUCT CODES

GM-08S01 Disposable Metzenbaum Monopolar Scissor 5mm X 330mm

GM-08S02 Disposable Metzenbaum Monopolar Scissor 5mm X 330mm

GM-08D01 Disposable Maryland Dissector 5mm X 330mm

GM-08G01 Disposable Lap Clinch with Ratchet 5mm X 330mm

GM-08G02 Fenestrated Grasper with Ratchet 5mm X 330mm

GM-08G04 Disposable Fenestrated Grasper (Long) with Ratchet 5mm X 330mm

GM-08G05 Disposable Babcock Grasper 5mm X 330mm

Disposable Clip Applier



DETAILS

Finger Rotation Wheel 360°

Tip Angled 20° To Facilitate Visualization Of Anatomic Features

20 Titanium Medium/Large Ligation Clips Are Pre-Loaded In The Device

29cm X 10mm Lucent Shaft Provides Line Of Sight To Remaining Clips

Ergonomic Hand Piece Facilitates Control

PRODUCT CODES

GM-0702 Disposable Clip Appliers 10mm x 330mm

Veress Needle



DETAILS

visual needle position indicator.

Sharp tip for easy insertion

High flow stopcock

Spring-loaded stylet provides audible "click" when passing through tissue layers.

PRODUCT CODES

GN-VN120 Veress needle 120mm (Audible click & indicator)

GM-VN150 Veress needle 150mm (Audible click & indicator)

Retrieval Bags



Endoscopic Retrieval Bag

DETAILS

Stainless steel ring open the bag for easy use.

High strength nylon material bag, prevent tissue fluid leakage and tears.

Easy to use, shorten the duration of operation

PRODUCT CODES

GM-0602 Disposable Endoscopic Retrieval Bags 100 x 150mm

GM-0603 Disposable Endoscopic Retrieval Bags 130mm x 180mm

GM-0604 Disposable Endoscopic Retrieval Bags 170mm x 240mm

GM-0605 Disposable Endoscopic Retrieval Bags 200mm x 260mm

Standard Retrieval Bag

DETAILS

Ripstop Nylon Material with Polymer Coating to Prevent Leakage

Nylon String for Ease & Comfort of Extraction

Various Sizes Available

PRODUCT CODES

GM-05BS Retrieval Bag Standard 80mm x 160mm

GM-05BM Retrieval Bags Standard 100mm x 150mm

GM-05BL Retrieval Bags Standard 170mm x 240mm



Suction Irrigation



DETAILS

Soft co-joined, low torque tubing - Minimizes intertwining for surgeons to work more freely.

Tubing is easy peel for separation.

Multi distal lateral holes to prevent blockage and help to irrigate every corner of abdominal cavity

Colour coded buttons to identify irrigation and suction.

Ergonomic handle for comfort. Surgeons feeling less fatigued especially for long and complexed procedures.

PRODUCT CODES

- GM-0421 Suction Irrigation 10mm x 330mm Shaft only
- GM-0422 Suction Irrigation Dual Spike 5mm x 330mm Shaft
- GM-0423 Suction Irrigation 5mm x 450mm Shaft only



GENMED
ENTERPRISES

CONTACT

For more information about our products and pricing,
please contact us:

✉ info@genmedltd.com

☎ 02920 361444

Quality is not a problem,
but a solution.