

Filters - For Buffalo Filter® Evacuators



VS/VET01



SING VSEY



VSEY1-4



PRFIL20

ViroSafe® Veterinary Filter

Designed to fit PlumeSafe® Turbo; Porta PlumeSafe® 604
Lasts for 24 hours of smoke evacuation time and captures 0.1 to 0.2 micron particles at 99.999% efficiency.

Dual male ports — 7/8 in (22 mm) and 1-1/4 in (32 mm)

Most economical choice for animal hospitals and clinics.

(1/case).....VSVET01

ViroSafe® Erbium: YAG Filter

Designed to fit PlumeSafe® Turbo; Porta PlumeSafe® 604
Specifically designed for use with the Erbium: Yag laser applications.
Lasts for 12 hours of smoke evacuation time when used with pre-filters BFEY04PF and captures particles 0.1 to 0.2 micron particles at 99.999% efficiency.

Dual male ports — 7/8 in (22 mm) and 1-1/4 in (32 mm).

(1/case).....SING VSEY

ViroSafe® Erbium: YAG Kit

Designed to fit Erbium: Yag laser applications

Using Erbium: YAG pre-filters helps to protect the efficacy of the ViroSafe® Erbium: YAG main filter (ensuring that it lasts for 12 hours of smoke evacuation time).

The pre-filter has a 1-3/8 in (35 mm) male port at each end—one to attach to the 1-1/4 in (32 mm) opening on the Main Erbium: YAG filter and another to attach to a 1-3/8 in (35 mm) tubing set.

Efficiency:

Main Filter: captures particles 0.1 to 0.2 micron particles at 99.999% efficiency

Pre-filter: 0.3–0.5 micron filtration

Usage Time

Main Filter: 12 hours

Pre-filter: 1–3 procedures

(1 Main Filter, 4 Pre-filters) VSEY1-4

Pre-filter

Not available in the United States

Pre-filter, Single Use

(20/case)..... PRFIL20