




Circular Stapler



Tilt top anvil for safe and easy removal through the anastomosis



Variable closed staple height to accommodate different tissue thickness



The cutting washer provides for positive audible and tactile feedback when fired



Inside and outside safety mechanism to prevent misfiring

SIZES AVAILABLE -

21mm | 24mm | 26mm | 29mm | 32mm

Product Code	Description
FCSME21	Circular stapler with tilt top, 21mm, dual safety, 4.5mm staple height
FCSME24	Circular stapler with tilt top, 24mm, dual safety, 4.5mm staple height
FCSME26	Circular stapler with tilt top, 26mm, dual safety, 4.8mm staple height
FCSME29	Circular stapler with tilt top, 29mm, dual safety, 4.8mm staple height
FCSME32	Circular stapler with tilt top, 32mm, dual safety, 5.0mm staple height



FERDINAND MEDICAL
Advanced Equipment, Trusted Care

QUALITY CARE THROUGH INNOVATION



Our Stapling Portfolio

Circular Stapler



Linear Cutter Stapler



Open Linear Stapler



Endo Stapling System

