



# FERDINAND MEDICAL

*Advanced Equipment, Trusted Care*

## LaparoVue® Visibility System Laparoscopic Solutions

Understanding a surgeon's need to see clearly during laparoscopic surgery, the LaparoVue Visibility System is an all-in-one solution designed to take the complexity out of achieving optimal visualization. An easy to use, disposable device that warms, white balances, cleans and defogs laparoscopes.



**AIRSEAL® iFS**  
ACCESS MANAGEMENT SYSTEM



**PLUMEPEN®**  
SMOKE EVACUATION PENCILS



**VISICLEAR®**  
SMOKE EVACUATOR



To learn more about these and other innovative products, visit:

[FERDINANDMEDICAL.COM](http://FERDINANDMEDICAL.COM)



# LaparoVue®

All-in-One Solution for Improved Laparoscopic Visibility

## Warms and White Balances


- Warms a laparoscope in 5 minutes to prevent fogging
- Scope Cradle supports and protects during warming and white balancing
- Designed for stability, LaparoVue eliminates surfactant leaks while warming and white balancing

## Cleans and Defogs

- Provides crisp, clear view in 5 seconds
- Offers a chemical and mechanical clean for superior residue removal while helping to prevent intra-operative laparoscope fogging
- Cleaning and defogging port designed to retain surfactant
- Radiopaque VueTip™ Trocar Swabs remove debris and excess liquid from trocar cannula
- Microfiber Cloth captures large residue and excess surfactant



# Distinctive By Design



## LaparoVue Product Features

### One Unit, Two Ports

- Warm a laparoscope while cleaning another, reduces preparation time when changing laparoscopes
- Accepts 3 – 12mm laparoscopes without the need for a port size converter
- Warms, white balances, cleans and defogs for up to 5 hours

### Designed for Stability

- Scope Cradle safely and securely supports laparoscope to reduce the risk of damage
- Upright functionality reduces chances of surfactant spillage

### Easy to remove batteries for disposal

### No assembly requirement streamlines procedure set-up time

### Single use for maximum infection control



**FERDINAND MEDICAL**  
*Advanced Equipment, Trusted Care*

**LAPAROVUE®  
VISIBILITY  
SYSTEM**

## Ordering Information

Description	Qty	Catalog Number
LaparoVue® Visibility System	10/Box	LAPVUE10
1 LaparoVue® Unit (Main Unit and Scope Cradle)		
2 VueTip™ Trocar Swabs (one large and one small)		
1 Microfiber Cloth		
VueTip™ Trocar Swab Accessory Kit	10/Box	VUETIP10
2 VueTip™ Trocar Swabs (one large and one small)		
1 Microfiber Cloth		