



**BUFFALO FILTER®**

# PlumePen® Pro

## Surgical Smoke Evacuation Pencil

The PlumePen® Pro features an easy to remove x-ray detectable capture port that allows for effortless electrode replacement to accommodate surgeries where longer electrodes are desired. PlumePen® Pro is smaller than competitive surgical smoke evacuation pencils and similar in size to standard electrosurgical pencils.



**SYSTEM 5000™**  
GENERATOR



**PLUMEPEN® ELITE**  
PENCILS



**PLUMEPEN® ULTRA**  
PENCILS



**VISICLEAR®**  
SURGICAL SMOKE EVACUATOR



To learn more about these and other innovative products, call

[FERDINANDMEDICAL.COM](http://FERDINANDMEDICAL.COM)

**BUFFALO FILTER®**

**Features & Benefits**

**Sleek, Slim Profile**

- Smaller than competitive surgical smoke evacuation pencils; similar in size to standard ESU pencils

**Removable Smoke Capture Port**

- Allows for easy electrode and capture port replacement

**Compatible with Standard Electrodes**

- Designed for surgeries where longer electrodes are desired

**Full Blade Tip Reveal**

- Provides an unobstructed view of surgical site

**360 Degree Swivel**

- Allows for fluid hand movement

**Textured Profile**

- Improves surgeon's grip



To learn more about this and other innovative products,  
[FERDINANDMEDICAL.COM](http://FERDINANDMEDICAL.COM)

Ordering Information

Description	Qty	Catalog Number
PlumePen® Pro (Non-Stick ESU Blade) 10 ft (3 m) tubing	20/Case	PLPR04020
PlumePen® Pro (Non-Stick ESU Blade) 15 ft (4.6 m) tubing	20/Case	PLPR04520
Capture Ports (X-ray detectable) Accommodates 4 in (102 mm) electrodes	20/Case	PROTIP420
Capture Ports (X-ray detectable) Accommodates 5 in (127 mm) electrodes	20/Case	PROTIP520
Capture Ports (X-ray detectable) Accommodates 6 in (152 mm) electrodes	20/Case	PROTIP620