



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

NEW

PlumePen® Ultra Surgical Smoke Evacuation Pencil

No matter the surgical procedure, PlumePen® Ultra is up to the challenge. An addition to the PlumePen portfolio, Ultra offers a telescopic capture port with an advanced swivel design that will delight the surgeon with superior hand movement and less wrist drag.



SYSTEM 5000™
GENERATOR



PLUMEPEN® ELITE
PENCILS



GOLDVAC® SLIM
PENCILS



VISICLEAR®
SURGICAL SMOKE EVACUATOR



To learn more about these and other innovative products,

FERDINANDMEDICAL.COM



Features & Benefits

Telescopic Capture Port with Enclosed Blade

- Blade moves in unison with capture port affording versatility in any surgical situation. Integrated design eliminates need for blade replacement saving time and money

Cam Lock System


- Prevents capture port movement during surgery. Audible click upon activation/deactivation of locking mechanism

Responsive, Multi-Directional Swivel

- Unique, floating ball joint design allows unrestricted hand movement

Integrated Wire/Tubing Management

- Ultra-light, PVC free, flexible tubing designed to reduce wrist pull and hand fatigue



Telescopic Capture Port with Enclosed Blade

To learn more about this and other innovative products,
FERDINANDMEDICAL.COM

Ordering Information

	Description	Qty	Catalog Number
NEW	PlumePen® Ultra (Non-Stick ESU Blade) with 3/8 in (9.5 mm) x 10 ft (3 m) hose/ tubing & 7/8 in (22mm) filter port adapter	20/Case	PLPUL2020
	PlumePen® Ultra (Non-Stick ESU Blade) with 3/8 in (9.5 mm) x 15 ft (4.6 m) hose/ tubing & 7/8 in (22 mm) filter port adapter	20/Case	PLPUL2520