

BUFFALO FILTER

PlumePort® SEO

Laparoscopic Smoke Filtration Device

PlumePort® SEO easily connects from a standard trocar to a smoke plume evacuator for active, efficient, and safe removal of surgical smoke during laparoscopic procedures. Utilizing available surgical smoke evacuators and coupling it with this flow controlled device, PlumePort® SEO provides fast and safe surgical smoke management while maintaining peritoneal distention without OR personnel interaction.



VISICLEAR®
SMOKE EVACUATION SYSTEM



LAPAROLIGHT®
VERESS NEEDLE



PLUMEPORT®
LAPAROSCOPIC



To learn more about these and other innovative products,

FERDINANDMEDICAL.COM

BUFFALO FILTER

Features & Benefits

Controlled Management

- Manages smoke throughout the surgical procedure

Universal Connection

- Standard luer lock fitting; connects to most surgical smoke evacuators

Up to 50% More Flow

- Provides greater flow than passive laparoscopic devices

Adjustable Flow Rate

- Allows flow adjustment via surgical smoke evacuator

Excellent Value

- An economical effective solution

Quiet and Easy to Use



To learn more about this and other innovative products,
FERDINANDMEDICAL.COM

Ordering Information

	Description	Qty	Catalog Number
NEW	PlumePort® SEO with 3/8 in (9.5 mm) x 10 ft (3 m) hose/ tubing & 7/8 in (22mm) filter port adapter	12/Case	PP012SEO