



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

Disposable Laparoscopic Specimen Bags

Self-opening, high-strength bags ensure safe specimen capture with easy, one-handed operation.



UNIVERSAL PLUS™
SUCTION/IRRIGATION



DETACHATIP® III
SYSTEM



VCARE®
UTERINE MANIPULATOR



To learn more about these and other innovative products, visit

FERDINANDMEDICAL.COM



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

DISPOSABLE
LAPAROSCOPIC
SPECIMEN BAGS

Sterile, disposable device designed for use as a receptacle for the collection and extraction of tissue, organs, and calculi during laparoscopic surgical procedures. The Specimen Bag System consists of a flexible memory wire specimen bag and introducer. The Specimen Bag System can be used through a 5, 10, or 12mm cannula.

Unique Design

The unique patented design ensures ease of use, and can be deployed and manipulated with one hand. The self-opening memory wire, high-strength large bag offers easy and safe specimen capture while resisting tearing, reducing the potential for cell migration.

Easy to Use

Indicators on the thumb ring position the opening (top) of the bag for correct placement. Bag can be detached within abdominal cavity for later retrieval.

Benefits

- Large bag volume to capture large specimens
- Minimum insertion force enables smooth deployment
- Ergonomic design provides easy one-handed operation
- Insert indicator ensures bag will be in proper position when deployed



Ordering Information

Description	Cannula	Diameter	Length	Volume	Qty Per Box	Catalog Number
Attached Memory Wire Bag, Polymer	5mm	3in (7.62cm)	4in (10.2cm)	120ml	10/case 1/pkg	SB534
Independent Memory Wire Bag, Polymer	10mm	3in (7.62cm)	6in (15.24cm)	255ml	10/case 1/pkg	SB836
Independent Memory Wire Bag, Polymer	10mm	5in (12.7cm)	7in (17.78cm)	680ml	10/case 1/pkg	SB857
Attached Memory Wire Bag, Polymer	10mm	3in (7.62cm)	6in (15.24cm)	255ml	5/case 1/pkg	SB936
Attached Memory Wire Bag, Polymer	10mm	5in (12.7cm)	7in (17.78cm)	680ml	5/case 1/pkg	SB957
Attached Memory Wire Bag, Polymer	12mm	7.5in (19.1cm)	9in (22.9cm)	1500ml	5/case 1/pkg	SB979
Attached Memory Wire Bag, Nylon	10mm	3in (7.62cm)	6in (15.24cm)	255ml	5/case 1/pkg	SB1036
Attached Memory Wire Bag, Nylon	10mm	5in (12.7cm)	7in (17.78cm)	680ml	5/case 1/pkg	SB1057
Attached Memory Wire Bag, Nylon	12mm	7.5in (19.1cm)	9in (22.9cm)	1500ml	5/case 1/pkg	SB1079