



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

Trocar


PRODUCT RANGE

Contents



- 1 Optical Trocars
 - 2 Bladeless Trocars
 - 3 Bladed Trocars
 - 4 Balloon Trocars
-

Optical Trocars



DETAILS

- Optical trocar design for safe insertion under vision
- Vision is maintained for the duration of insertion
- Patent protected seal technology
- Patent protected lock mechanism to fix scope securely
- Dilating wings to assist with ease of insertion whilst maintaining maximum control.
- Bladeless tip to minimise port site injury and reduce rate of port site herniation
- Threaded and balloon cannula available

PRODUCT CODES

GN-TO0575	Bladeless Threaded Optical Trocar 5mm x 75mm
GN-TO0510	Bladeless Threaded Optical Trocar 5mm x 100mm
GN-TO0515	Bladeless Threaded Optical Trocar 5mm x 150mm
GN-TO1010	Bladeless Threaded Optical Trocar 10mm x 100mm
GN-TO1210	Bladeless Threaded Optical Trocar 12mm x 100mm
GN-TO1215	Bladeless Threaded Optical Trocar 12mm x 150mm
GN-TO1510	Bladeless Threaded Optical Trocar 15mm x 100mm

* Dual packs and individual cannulas also available

Bladeless Trocars



DETAILS

- Ergonomic design for safe insertion
- Threaded cannula for stability of trocar
- Patent protected seal technology
- Dilating wings to assist with ease of insertion whilst maintaining maximum control.
- Bladeless tip to minimise port site injury and reduce rate of port site herniation
- Threaded and balloon cannula available

PRODUCT CODES

GN-TW0575	Bladeless Threaded Trocar 5mm x 75mm
GN-TW0510	Bladeless Threaded Trocar 5mm x 100mm
GN-TW0515	Bladeless Threaded Trocar 5mm x 150mm
GN-TW1010	Bladeless Threaded Trocar 10mm x 100mm
GN-TW1210	Bladeless Threaded Trocar 12mm x 100mm
GN-TW1215	Bladeless Threaded Trocar 12mm x 150mm
GN-TW1510	Bladeless Threaded Trocar 15mm x 100mm

* Dual packs and individual cannulas also available

Bladed Trocars



DETAILS

- Ergonomic design for safe and accurate insertion
- Integrated shield reset button
- Blade shield protects internal structures from damage during insertion
- Patent protected seal technology
- High grade steel blades and rods are designed to assist with dilating as well as cutting •

Low level insertion force required in turn preventing damage to internal structures

- Threaded and balloon cannula available

PRODUCT CODES

GN-TB0575	Bladed Threaded Trocar 5mm x 75mm
GN-TB0510	Bladed Threaded Trocar 5mm x 100mm
GN-TB1010	Bladed Threaded Trocar 10mm x 100mm
GN-TB1210	Bladed Threaded Trocar 12mm x 100mm

* Dual packs and individual cannulas also available

Balloon Trocars



DETAILS

- Various style obturators available – Optical, Bladeless, Blunt, Bladed
- Unique fixation device for stronger stability
- Patent protected seal technology
- Patent protected indicator balloon for gauging pressure in primary balloon
- Smooth cannula
- Ring fixation and Anchor fixation devices available

PRODUCT CODES (ANCHOR FIXATION DEVICE)

- | | |
|-------------|--|
| GN-TOHF0510 | Blunt Balloon Trocar 5mm x 100mm with anchor fixation |
| GN-TOHF1010 | Blunt Balloon Trocar 10mm x 100mm with anchor fixation |
| GN-TOHF1210 | Blunt Balloon Trocar 12mm x 100mm with anchor fixation |

PRODUCT CODES (ANCHOR FIXATION DEVICE)

- | | |
|-----------|--|
| GN-TH0510 | Blunt Trocar 5mm x 100mm with anchor fixation |
| GN-TH1010 | Blunt Trocar 10mm x 100mm with anchor fixation |
| GN-TH1210 | Blunt Trocar 12mm x 100mm with anchor fixation |
-

Balloon Trocars



PRODUCT CODES (RING-FIXATION DEVICE)

GN-TWF0575	Bladeless Balloon Trocar 5mm x 75mm with ring fixation
GN-TWF0510	Bladeless Balloon Trocar 5mm x 100mm with ring fixation
GN-TWF0515	Bladeless Balloon Trocar 5mm x 150mm with ring fixation
GN-TWF1010	Bladeless Balloon Trocar 10mm x 100mm with ring fixation
GN-TWF1210	Bladeless Balloon Trocar 12mm x 100mm with ring fixation
GN-TWF1215	Bladeless Balloon Trocar 12mm x 150mm with ring fixation
GN-TWF1510	Bladeless Balloon Trocar 15mm x 100mm with ring fixation

PRODUCT CODES (RING-FIXATION DEVICE)

GN-TOF0575	Optical Balloon Trocar 5mm x 75mm with ring fixation
GN-TOF0510	Optical Balloon Trocar 5mm x 100mm with ring fixation
GN-TOF0515	Optical Balloon Trocar 5mm x 150mm with ring fixation
GN-TOF1010	Optical Balloon Trocar 10mm x 100mm with ring fixation
GN-TOF1210	Optical Balloon Trocar 12mm x 100mm with ring fixation
GN-TOF1215	Optical Balloon Trocar 12mm x 150mm with ring fixation
GN-TOF1510	Optical Balloon Trocar 15mm x 100mm with ring fixation

PRODUCT CODES (RING-FIXATION DEVICE)

GN-TBF0575	Bladed Balloon Trocar 5mm x 75mm with ring fixation
GN-TBF0510	Bladed Balloon Trocar 5mm x 100mm with ring fixation
GN-TBF1010	Bladed Balloon Trocar 10mm x 100mm with ring fixation
GN-TBF1210	Bladed Balloon Trocar 12mm x 100mm with ring fixation



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

Quality is not a problem,
but a solution.