

BUFFALO FILTER

ViroVac® Surgical Smoke Evacuator


ViroVac® surgical smoke plume evacuator provides ideal flow for physician offices and outpatient surgery centers. It combines customizable (variable) filter life, a compact shelf design, 3-port filter with filter life tracking technology, while exceeding expectations for value and efficacy.



PLUMEPEN® ELITE
PENCILS



PLUMEPEN® ULTRA
PENCILS



GOLDVAC® SLIM
PENCILS



To learn more about these
and other innovative
products, visit
FERDINANDMEDICAL.COM

BUFFALO FILTER

VIROVAC®
SURGICAL SMOKE
EVACUATOR

Features & Benefits

Compact Design

- Allows for easy placement on shelf or counter and accommodates stackability.

Patented 3-Port Filter Design


- Versatile 3-port filter design features covered, hidden male ports and accommodates a variety of tubing sizes providing versatility in procedure and accessory setup.

Variable Filter Life (up to 35 hours)

- Delivers unsurpassed value and reduces per procedure costs.

Filter Tracking Technology

- Automatically tracks and records filter life through a patented technology.



To learn more about this and other innovative products, visit FERDINANDMEDICAL.COM

Ordering Information

Description	Qty	Catalog Number
ViroVac® (100/120 Volt Systems) Includes one filter and one footswitch	1/Each	W120
ViroVac® (200/240 Volt Systems) Includes one filter and one footswitch	1/Each	VV220
EZLink® Sold separately	1/Each	EZLINK01
ViroSafe® Fluid Trap	10/Case	VSFT10
ViroSafe® ViroVac® Filter	1/Case	VS353
	2/Case	VS35302
	4/Case	VS35304