



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

QUALITY **CARE** THROUGH **INNOVATION**



Product Catalog



Panther Healthcare USA



PANTHER HEALTHCARE has operated for 22 years and is the 3rd largest global manufacturer of open and endoscopic staplers. We distribute our portfolio in over 80 countries. Panther Healthcare USA is proud to bring our innovative portfolio of products to the United States.

Every surgical stapler from Panther Healthcare is subjected to a firing test before being packaged and sterilized for patient use. Every single stapler and cartridge is photographed and stored for QMS tracking.

OUR VISION:

The Company that always exceeds your expectations.

OUR VISION:

- To provide high quality products with an innovative future that connects value and fiscal responsibility to the healthcare community.
- We will achieve our mission by:
 - Continuously focusing on product/process improvement and innovation in all our offerings
 - Providing the most cost-effective solution on every contract
 - Relentless pursuit of best-in-class order fulfillment

VALUE PROPOSITION:

Innovation, Quality and Consistency at a Great Price.

Contents

STAPLING PORTFOLIO

Endo Staplers

- 2 Endo Linear Cutter Handles
- 3 Endo Linear Cutter Cartridges

Open Staplers

- 4 Linear Cutter and Cartridges
- 5 Linear Stapler and Cartridges
- 6 Circular Stapler
- 6 Extra Long Circular Stapler
- 7 Linear Stapler with Knife and Cartridges
- 7 Panther Skin Stapler & Extractor

- 8 Standard Terms & Conditions
- 9 Contact Information

Endo Linear Cutter Handles

The **PANTHER Manual Endo Linear Cutter Handle** is available in 3 different shaft lengths, Short, Standard and XL. The Panther Endo Linear Cutter simultaneously staples and divides tissue while placing two triple staggered rows of titanium staples. The universal stapler handle is designed to accommodate all cartridges and can be fired up to 25 times per patient.

ENDO MANUAL LINEAR CUTTER HANDLE



Product Code	Description	Shaft Length	Pack Qty.
CEAC 30	Endo Linear Cutter (Standard Handle)	160mm	6 / box
CEAC 45	Endo Linear Cutter (Short Handle)	60mm	6 / box
CEAC 60	Endo Linear Cutter (Extra Long Handle)	260mm	6 / box

The **PANTHER SMART Endo Linear Cutter Handle** is a battery powered device that includes an RFID reader that alerts the stapler to changing tissue thickness. The SMART Handle is available in 3 different shaft lengths-Short, Standard and XL. The Endo Linear Cutter Stapler simultaneously staples and divides tissue while placing two triple staggered rows of titanium staples. The SMART handle is designed to accommodate all cartridges and can be fired up to 15 times per patient.

SMART ENDO POWERED LINEAR CUTTER HANDLE



Product Code	Description	Shaft Length	Pack Qty.
PEAS	Panther Smart Handle (Short Handle)	60mm	6 / box
PEAM	Panther Smart handle (Standard Handle)	160mm	6 / box
PEAL	Panther Smart Handle (Extra Long Handle)	260mm	6 / box

Endo Linear Cutter Cartridges

P2G ARTICULATING CARTRIDGE WITH DOUBLE GRIP TECHNOLOGY



	Product Code	Staple-Length	Staple-Height	Pack Qty.
○	CADF-30T	30mm	2.5mm	6 / box
●	CADF-30N	30mm	3.5mm	6 / box
○	CADF-45T	45mm	2.5mm	6 / box
●	CADF-45N	45mm	3.5mm	6 / box
●	CADF-45S	45mm	4.2mm	6 / box
●	CADF-45D	45mm	4.8mm	6 / box
●	CADF-45R	45mm	5.0mm	6 / box
○	CADF-60T	60mm	2.5mm	6 / box
●	CADF-60N	60mm	3.5mm	6 / box
●	CADF-60S	60mm	4.2mm	6 / box
●	CADF-60D	60mm	4.8mm	6 / box
●	CADF-60R	60mm	5.0mm	6 / box

P2G CURVED TIP ARTICULATING CARTRIDGE WITH DOUBLE GRIP TECHNOLOGY



	Product Code	Staple-Length	Staple-Height	Pack Qty.
○	CADJ-30T	30mm	2.5mm	6 / box
●	CADJ-30N	30mm	3.5mm	6 / box
○	CADJ-45T	45mm	2.5mm	6 / box
●	CADJ-45N	45mm	3.5mm	6 / box
●	CADJ-45S	45mm	4.2mm	6 / box
●	CADJ-45D	45mm	4.8mm	6 / box
●	CADJ-45R	45mm	5.0mm	6 / box
○	CADJ-60T	60mm	2.5mm	6 / box
●	CADJ-60N	60mm	3.5mm	6 / box
●	CADJ-60S	60mm	4.2mm	6 / box
●	CADJ-60D	60mm	4.8mm	6 / box
●	CADJ-60R	60mm	5.0mm	6 / box

Linear Cutter Stapler with Extra Grip Technology

The **Panther Linear Cutter Stapler** places two double staggered rows of titanium staples while simultaneously dividing the tissue between the 2 rows of staples. The staplers and reloads are available in 2 different lengths - 60mm and 80mm and 2 staple sizes - 3.8mm (blue) and 4.8mm (Green). The stapler can be reloaded up to 7 times for a total of 8 firings in a single procedure

LINEAR CUTTER STAPLER



	Product Code	Description	Staple-Length Height		Pack Qty.
●	SSAB-60 3.8MM	Linear Cutter Stapler	60mm	3.8mm	6 / box
●	SSAB-60 4.8mm	Linear Cutter Stapler	60mm	4.8mm	6 / box
●	SSAB-80 3.8mm	Linear Cutter Stapler	80mm	3.8mm	6 / box
●	SSAB-80 4.8mm	Linear Cutter Stapler	80mm	4.8mm	6 / box

CARTRIDGES FOR LINEAR CUTTER STAPLER



	Product Code	Staple-Length	Staple-Height	Pack Qty.
●	SADB-60N	60mm	3.8mm	6 / box
●	SADB-60D	60mm	4.8mm	6 / box
●	SADB-80N	80mm	3.8mm	6 / box
●	SADB-80D	80mm	4.8mm	6 / box

Linear Stapler

The **Panther Linear Stapler** places a double staggered row of titanium staples. Our exclusive color identification system located at the back of the handle, identifies each phase of the firing cycle. The stapler reloads are available in lengths of 30mm, 45mm, 60mm and 90mm with 2 different staple heights - 3.5mm (blue) and 4.8mm (green). The stapler can be reloaded up to 7 times for a total of 8 firings in a single procedure.

LINEAR STAPLER



	Product Code	Description	Staple-Length Height		Pack Qty.
●	FLSLF30-3.5	Linear Stapler	30mm	3.5mm	3 / box
●	FLSLF30-4.8	Linear Stapler	30mm	4.8mm	3 / box
●	FLSLF45-3.5	Linear Stapler	45mm	3.5mm	3/ box
●	FLSLF45-4.8	Linear Stapler	45mm	4.8mm	3 / box
●	FLSLF60-3.5	Linear Stapler	60mm	3.5mm	3/ box
●	FLSLF60-4.8	Linear Stapler	60mm	4.8mm	3 / box
●	FLSLF90-3.5	Linear Stapler	90mm	3.5mm	3/ box
●	FLSLF90-4.8	Linear Stapler	90mm	4.8mm	3 / box

CARTRIDGES FOR LINEAR STAPLER



	Product Code	Staple-Length	Staple-Height	Pack Qty.
●	FLSLH-3.5A	30mm	3.5mm	3 / box
●	FLSLH-4.8A	30mm	4.8mm	3 / box
●	FLSLH-3.5B	45mm	3.5mm	3/ box
●	FLSLH-4.8B	45mm	4.8mm	3 / box
●	FLSLH-3.5C	60mm	3.5mm	3/ box
●	FLSLH-4.8C	60mm	4.8mm	3 / box
●	FLSLH-3.5E	90mm	3.5mm	3/ box
●	FLSLH-4.8E	90mm	4.8mm	3 / box

Circular Stapler & Extra Long Circular Stapler

The Panther Circular Stapler places a circular, double staggered row of titanium staples while resecting the tissue and creating a circular anastomosis. Our exclusive combination of variable staple height closure and tilt top anvil delivers a secure anastomosis over a wide range of tissue, with easy and safe removal post firing. The Circular Staplers are available in 5 different diameters – 21mm, 24mm, 26mm, 29mm and 32mm.



CIRCULAR STAPLER

	Product Code	Description	Staple Height	Closed Staple Height Range	Pack Qty.
●	FCSME21	Circular Stapler with Tilt Top Anvil, 21mm	4.5mm	1.4-2.2mm	3 / box
●	FCSME24	Circular Stapler with Tilt Top Anvil, 24mm	4.5mm	1.4-2.2mm	3 / box
○	FCSME26	Circular Stapler with Tilt Top Anvil, 26mm	4.8mm	1.6-2.4mm	3 / box
●	FCSME29	Circular Stapler with Tilt Top Anvil, 29mm	4.8mm	1.6-2.4mm	3 / box
●	FCSME32	Circular Stapler with Tilt Top Anvil, 32mm	5.0mm	1.8-2.6mm	3 / box

The Panther Disposable Extra Long Circular Staplers place a circular double staggered row of titanium staples while resecting the tissue and creating a circular anastomosis. The stapler is hermetically sealed for laparoscopic procedures and is available in 25mm, 29mm and 32mm diameters.



EXTRA LONG CIRCULAR STAPLER

	Product Code	Description	Diameter	Staple Height	Closed Staple Height Range	Pack Qty.
	FCSLWAF25	Extra Long Circular Stapler w/ Tilt Top Anvil	25mm	4.8mm	1.6-2.4mm	3 / box
	FCSLWAF29	Extra Long Circular Stapler w/ Tilt Top Anvil	29mm	4.8mm	1.6-2.4mm	3 / box
	FCSLWAF32	Extra Long Circular Stapler w/ Tilt Top Anvil	32mm	5.0mm	1.8-2.6mm	3 / box

Linear Stapler with Knife

The Panther Linear Stapler with Knife places two double rows of staples and cuts in between them. The cutting direction is perpendicular to the stapling one: first stapling then cutting ensures a smoother cut on thicker tissue.



CAZB-40D
Reloading Unit

LINEAR STAPLER WITH KNIFE

	Product Code	Description	Staple Line Length	Cutting Line Length	Staple Height	Staple Formation	Pack Qty.
	CAAB-40	Disposable Linear Cutter Stapler	43mm	40mm	4.8mm	2.0mm	3 / box
	CAZB-40D	Reloading Unit for Linear Cutter Stapler	43mm	40mm	4.8mm	2.0mm	1 / box

Single-Patient Use Skin Stapler & Extractor

The Panther Single-Patient Use Skin Stapler is designed to be fired from multiple angles. It has an ejector spring that automatically disengages the staple to facilitate instrument removal in any direction.



DSEB
Staple Extractor

SINGLE-PATIENT USE SKIN STAPLER

	Product Code	Description	Wire Diameter	Staple Width	Staple Height	Pack Qty.
	DSSA-35W	Single-Patient Use Skin Stapler	.6mm	12.8mm	3.5mm	6 / box
	DSEB	Single-Patient Use Staple Extractor	—	—	—	3 / box

For product use and safety information, please refer to the instructions for use manual.

Standard Terms & Conditions

ORDER PLACEMENT:

All purchase orders must be emailed to orders@pantherhealthcareusa.com or faxed to 646-858-0166. All purchase orders must be in typed format, sent as a PDF file, on purchasing entity letterhead.

ORDER SHIPPING:

All purchase orders received before 1:00pm EST will be shipped out same day by Standard ground. At time of shipping, all orders will be charged a shipping charge. Next day delivery rush orders will be charged the Expedited Shipping fee at time of shipping.

PAYMENT TERMS:

30 days Net on invoice from date of invoice. No statement will be issued. 1.5% monthly interest on invoices not paid within 30 days of original invoice date.

PRODUCT RETURNS:

Please contact orders@pantherhealthcareusa.com to request a Panther Healthcare USA Return Goods Authorization Number (RGA#). Any product that is shipped to Panther Healthcare USA without a Panther RGA# issued will not be accepted. When an RGA# is obtained, please ship the product pre-paid to Panther Healthcare USA with the RGA# clearly marked on outer package(s) containing the product return. All products being returned with an RGA# must be in sterile, re-sellable condition, in original outer packaging, and a copy of the original invoice must be attached. Refunds for returned products are paid at Panther Healthcare USA's discretion. A 20% restocking fee will be charged if the product is returned within 30 days of original purchase date. No credit will be given if the product is returned after 30 days of the original purchase date. Products being returned because of claimed defects are not covered by this returns policy, but instead are subject to Panther Healthcare USA's limited warranty.

LIMITED WARRANTY:

Panther Healthcare USA ("Seller") warrants that the products sold hereunder will be free from defects in material and workmanship until the date of expiration provided on the product. Any claim made under this limited warranty must be received by Seller (productcomplaints@pantherhealthcareusa.com) within thirty (30) days of the product's date of expiration. SELLER'S WARRANTY HEREIN IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES OF SELLER, ITS PARENT AND AFFILIATED ENTITIES AND THE MANUFACTURER OF THE PRODUCTS (COLLECTIVELY, "Manufacturing and Selling Parties"), WHETHER EXPRESS, IMPLIED, STATUTORY, OR OTHERWISE CREATED UNDER APPLICABLE LAW, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY OF MERCHANTABILITY AND ANY WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE OR USE. IN NO EVENT SHALL SELLER OR THE MANUFACTURING AND SELLING PARTIES BE LIABLE (AND PURCHASER SHALL NOT ASSERT ANY CLAIM) FOR SPECIAL, INCIDENTAL, PUNITIVE, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING LOSS OF PROFITS. PURCHASER'S SOLE AND EXCLUSIVE REMEDIES UNDER THIS WARRANTY SHALL BE LIMITED, AT SELLER'S EXCLUSIVE DISCRETION, TO: (1) REPLACEMENT OF ANY DEFECTIVE PRODUCTS OR PART THEREOF; OR (2) REPAIR OF ANY DEFECTIVE PRODUCTS OR PART THEREOF; OR (3) RETURN OF THE PURCHASE PRICE FOR ANY DEFECTIVE PRODUCT OR PART THEREOF. The warranty and remedy set forth above are conditioned upon the proper storage, use and maintenance of the products and conformance with all applicable recommendations of the Manufacturing and Selling Parties with respect to the products. No agent, employee or representative of Seller (or any distributor, dealer or sales representative of Seller) has the authority to bind Seller to any affirmation, representation or warranty concerning the products sold hereunder, and unless such affirmation, representation or warranty is specifically included in these Terms and Conditions, it will not form a part of the basis of these Terms and Conditions and shall in no way be binding upon the Seller or enforceable by any purchaser of the products.

NOTES: This policy is subject to change without prior notice.

Panther Healthcare USA



CONTACT INFORMATION

FERDINANDMEDICAL.COM



**INNOVATION, QUALITY
& CONSISTENCY**
at a GREAT PRICE



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

FERDINANDMEDICAL.COM