



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

NEW

BUFFALO FILTER

ViroSafe® Fluid Trap Surgical Smoke Evacuation Accessories

Designed to capture fluids to help maintain the integrity of the filtration system and assist with maintaining filter life.



VISICLEAR®
SURGICAL SMOKE EVACUATOR



PLUMEPEN® ULTRA
PENCIL



PLUMEPEN® ELITE
PENCIL



SYSTEM 5000®
GENERATOR



To learn more about these and other innovative products,
FERDINANDMEDICAL.COM



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

VIROSAFE®
FLUID TRAP
SURGICAL SMOKE
EVACUATION
ACCESSORY

BUFFALO FILTER

Features & Benefits

Translucent, Spill Resistant Design

Simple to Read Measurement Increments

- (holds up to 100cc fluid capacity)

Easy to Attach and Remove

Connection Ports

- 7/8 in (22 mm) and 1/4 in (6.4 mm)

Fits Buffalo Filter ViroSafe® Filters

Cap Posts for Storage Convenience



To learn more about this and other innovative products, call **1-800-448-6506** or visit **CONMED.com**.

Ordering Information

	Description	Qty	Catalog Number
NEW	ViroSafe® Fluid Trap	10/Case	VSFT10

CONMED Corporation
525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-305

©2019 CONMED Corporation,
MC20190472 4/19