



Series	HC29	HO35	HC35	HO42
Insertion tube outer diameter	2.9 mm	3.5 mm	3.5 mm	4.2 mm
Instrument channel width	/	1.8 mm	/	
Compatibility	/	5 Fr	/	
Infusion	Single infusion	Continuous infusion		
Pixel	160 k		920 k	
Field of view	120 °		140 °	
Direction of view	15 °			
Depth of field	3-100 mm			
Working length	220 mm			
Weight	50 g			

Reference:
*Angelis, C. D. 2003. "Office Hysteroscopy and Compliance: Mini-Hysteroscopy versus Traditional Hysteroscopy in a Randomized Trial." Human Reproduction 18(11): 2441-45. doi:10.1093/humrep/deg463.

Healing for life, and life for healing

- Minimizing Environmental Impact
Our disposable endoscopes are made with environmental care.
- Caring for the Environment, Caring for You
Let's invest in a sustainable future together.



SeleneView


Single-Use Digital Hysteroscope ▼



H series

HC29 | HC35 | HO35 | HO42

Office | Outpatient | Ambulatory Surgery Room ▼



“See and Treat”
with single use hysteroscope



Lower Cost
Shorten hospital stay and reduce hospitalization costs for patients.




High Efficiency
Same-day surgery and discharge




Less Pain
Reduce cervical trauma and preserve fertility

Gentle Care ▼

Comfort in every moment




Vaginoscopic approach offers three key advantages:
No need for cervical dilation, no anesthesia, and no speculum.





Superior Visibility

Clarity in a new way ▼



Chromvue



HD HD Pixel with 140° Wide FOV

Optimized Efficiency

Convenient in every step ▼

Ultra light 01
Only 50 g weight, comfortable for whole day operation

See and treat 02
Compatible with 5 Fr instrument

Take and use 03
No need for assembling

5s 5s Fast uterine distension 04
Achieve fast clear vision