



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

www.ferdinandmedical.com



ENDO RLCZ

Endoscopic Linear Cutting Staplers



New
Design



Omnia Reload
Compatible




Improved
Force-to-fire



CE Certified

ENDO REACH RLCZ Universal Stapler

- Compatible with Omnia reloads
- Improved handle design



Product Code	Shaft Length (mm)	Max number of firing	QTY/ BOX
ENDO RLCSZ	85	25	3
ENDO RLCZ	155		
ENDO RLCLZ	255		

Compatible Reload Specifications

Product Code	Staple Line Length (mm)	Trocar Compatibility (mm)	Color	Closed Staple Height (mm)	Tissue Thickness	Number of staple rows	QTY/ BOX
ID30TAN	30	12	Tan	0.75/1.0/1.25	Vascular/ Thin	6	6
ID45TAN	45						
ID60TAN	60						
ID30PUL	30		Purple	1.25/1.5/1.75	Medium/ Thick		
ID45PUL	45						
ID60PUL	60						
ID45BLK	45	15	Black	1.75/2.0/2.2	Extra Thick		
ID60BLK	60						

Contacts

Ferdinand Medical ApS
8600 Silkeborg
Denmark

+45 80 29 10 00
info@ferdinandmedical.dk
www.ferdinandmedical.com

Importer



Distributor

