



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

www.ferdinandmedical.com



Single-Use Semi-Rigid Instrument for Hysteroscopy

Biopsy Forceps · Alligator Grasper · Endoscopic Scissors



1.6 mm / 5 Fr
Outer diameter



380 mm
Working length



115°
Jaw opening angle



Semi-Rigid
Single-Use

Single-Use Semi-Rigid Instrument for Hysteroscopy

Biopsy Forceps | Alligator Grasper | Endoscopic Scissors



Large jaw opening angle of 115° for sufficient tissue sampling




Unique inner jaw surface texture boosts clamping force



Smooth cutting with moderate opening and closing force

PISTOL-GRIP DESIGN

1.6 mm / 5 Fr outer diameter
380 mm working length
Semi-rigid shaft



Ordering information

Model	Description	Rigid	O.D.	Working length
BF-S-1.6x380	Single-Use Biopsy Forceps	Semi-Rigid	1.6 mm/5Fr	380 mm
GF-S-1.6x380	Single-Use Alligator Grasper	Semi-Rigid	1.6 mm/5Fr	380 mm
PS-S-1.6x380	Single-Use Endoscopic Scissors	Semi-Rigid	1.6 mm/5Fr	380 mm

Contacts

Ferdinand Medical ApS

8600 Silkeborg
Denmark

+45 80 29 10 00
info@ferdinandmedical.dk
www.ferdinandmedical.com