



FERDINAND MEDICAL

Advanced Equipment, Trusted Care

www.ferdinandmedical.com



IRIS Powered Stapler

Redefining the Reposable Stapling System



Tissue
Insight



Surgical
Reach



Optimized
Ergonomics



Reprocessing
Workflow





CE Certified

Tissue Insight

Visible Feedback Informed Decisions

Real Time Firing Curve™

From segmented indication to continuous insight




Real-time Audible & Visual Feedback

See. Hear. Act with Confidence.




Tissue condition alert

Prevents unsafe firing




Pre-compression

For optimal staples formation



Articulation Awareness

Alert at 0 & MAX angle



Surgical Reach


Designed for Access, Precision, and Adaptability.

1 Conventional Staplers

45°–55° Articulation Limited
Access in Confined Anatomy


2 Ultraflex 90

Up to 90° Articulation
Enabling Precise Horizontal
& Lower Distal Transection



Ultraflex™ 90

Extending Access Where it Matters Most



Curved Tip



Progressive Height



Optimized Ergonomics

Optimized Weight and Balance for Precise Control



Direction-True™ Articulation

Consistent Control at Any Orientation

One-Handed Control

No Re-mapping Required



Simplify Reprocessing Workflow



Competitor

Reusable

Handle build in Battery

- Charging downtime
- Multiple handles to coordinate

Re-sterilization



Reusable Adapter + Separate Shell

- Re-sterilization required
- Operational complexity

Disposable



IRIS

Reusable

Detachable Power System

- Quick battery swap
- One handle is enough

Fully Disposable

UniShell™ Adapter

- Single-Step Setup
- No re-sterilization

Powered by Sustainable Platform



IRIS Powered Stapler

A reusable powered stapling system combining Ultraflex™ 90 articulation, Real Time Firing Curve™ feedback and a detachable power system — designed for access, confidence and a faster workflow.



Disposable Instrument Product Codes

Product Code	Description	Instrument Length (mm)	Shaft Length (mm)
IDRL	Linear Stapler, Long	500	250
IDRM	Linear Stapler, Medium	400	150
IDRS	Linear Stapler, Short	350	100

Reusable Instrument Product Codes

Product Code	Description
IDRCORE	Handle
IDRCORE-0	Handle (without reusable battery pack)

Accessory Product Codes

Product Code	Description
IDRBATTERY	Reusable Battery Pack
IDRCHARGER	Battery Charger
IDRTOOL	Manual Override Tool

Contacts

Ferdinand Medical ApS
 8600 Silkeborg
 Denmark

 +45 80 29 10 00
 info@ferdinandmedical.dk
 www.ferdinandmedical.com

Importer



Distributor

