




FERDINAND MEDICAL

Advanced Equipment, Trusted Care

www.ferdinandmedical.com

- HIGH QUALITY
- PRECISION
- RELIABLE



Laparoscopic instruments for perfect performance

Monopolar Hand Instruments for laparoscopic surgery
Single use



5 mm
Port size



330 mm
Shaft length



Sterile
Single use



Box of 12
Packaging

Monopolar Laparoscopic Hand Instruments

Single use | For 5 mm ports | Shaft length 330 mm

The CIMPAX® Monopolar Laparoscopic Hand Instruments offer everything the professional physician requires to ensure perfect performance throughout the complete procedure.



Product features

- Lightweight well-balanced instruments
- Comfortable, easy open and close handle
- Naturally located lock for ratchet feature
- Free and simple rotation of jaws
- Double hinge design ensures easy and precise jaw motion and maximum opening
- For 5 mm ports. Shaft length 330 mm

HIGH QUALITY

PRECISION

RELIABLE



JAW DESIGN

Double hinge design ensures easy and precise jaw motion and maximum opening

Instrument range



C-CUT 11 or 17 mm blades



C-DISSECT 20 mm jaw length



C-GRASP 18 mm jaw length



C-GRASP 22 mm jaw length



C-CLINCH 22 mm jaw length



C-BAB 22 mm jaw length

Complete monopolar range

Scissors, dissector, atraumatic graspers, clinch grasper and Babcock forceps – all with 330 mm shaft for 5 mm ports, free rotation of jaws and a naturally located ratchet lock. Supplied sterile for single use, 12 pieces per box.

Metzenbaum

Maryland


Atraumatic

Clinch

Babcock

Ordering information

Product	Description	Size	Packaging	Ref.
C-CUT	Laparoscopic Scissors, Metzenbaum Curved	17 mm blades	Box of 12, sterile	C23-11-331
C-CUT	Laparoscopic Scissors, Metzenbaum Curved	11 mm blades	Box of 12, sterile	C23-11-337
C-DISSECT	Laparoscopic Dissector, Maryland Curved	20 mm jaws	Box of 12, sterile	C23-11-334
C-GRASP	Laparoscopic Grasper, Atraumatic with Ratchet	18 mm jaws	Box of 12, sterile	C23-11-332
C-GRASP	Laparoscopic Grasper, Atraumatic with Ratchet	22 mm jaws	Box of 12, sterile	C23-11-336
C-CLINCH	Laparoscopic Grasper, Clinch Grasper with Ratchet	22 mm jaws	Box of 12, sterile	C23-11-333
C-BAB	Laparoscopic Forceps, Babcock with Ratchet	22 mm jaws	Box of 12, sterile	C23-11-335



Monopolar Laparoscopic Hand Instruments

Single use · Sterile · For 5 mm ports · Shaft length 330 mm



Contacts

Ferdinand Medical ApS

8600 Silkeborg
Denmark

+45 80 29 10 00
info@ferdinandmedical.dk
www.ferdinandmedical.com